lumie[®] Desklamp[™]

OPERATING INSTRUCTIONS

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BRUKSANVISNING

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OPERATING INSTRUCTIONS

Desklamp

- 1 Introduction
- 2 Precautions
- 3 Safety
- 4 Getting started
- 5 Using Desklamp
- 6 When should I use it?
- 7 What time of day?
- 8 LED lifetime
- 9 Warranty
- 10 Technical specification
- 11 Contact us

1 Introduction

Using Desklamp will provide you with the bright light your body needs during the dark winter months. It can help to boost your energy levels, put you in a better mood and make you feel more awake¹. Light therapy is safe and is proven to help alleviate the winter blues and SAD (Seasonal Affective Disorder)².

Our circadian rhythms are strongly influenced by light and the short wavelengths around the blue part of the spectrum are known to be particularly important. Because the white LEDs used in Desklamp emit extra blue light, they offer an effective alternative to the white fluorescent tubes used in larger lightboxes.

Please read these instructions carefully to get the most out of your Desklamp.

1 J Affect Disord. 2000 Jan-Mar;57(1-3):55-61.

2 Am J Psychiatry. 2005 Apr;162(4):656-62.

	Serial number
	European Declaration of Conformity: Class IIa Medical Device
	Manufacturer: Lumie, 3 The Links, Bar Hill, Cambridge, CB23 8UD, United Kingdom.

	Do not place in household waste. Dispose of at a recycling centre.
	Read instructions for use.
	Caution

	Fuse
	Double-insulated (Class II electrical)
	AC (alternating current) power supply

	Member of the 'Green Dot' recycling scheme.
	Recyclable

2 Precautions

Research shows that light therapy is very safe. However, please check with your doctor first if any of the following apply:

- you have a major eye problem or have had eye surgery;
- you have been advised to avoid bright light or are taking medication that may cause photo-sensitivity;
- you have been diagnosed with depression or other mood/sleep disorders. This is especially important if you are also taking medication for this, since light therapy can also affect your mood.

A few people experience side effects but these are generally mild (including headaches, eye strain, nausea and hyperactivity) and can usually be solved simply by switching off the light. Over the next few days, gradually build up the time you spend in front of it to find a treatment time that works best for you. Some people complain of difficulty sleeping and often this is down to using the light too close to bedtime.

If you have any problems while using Desklamp, or if you've been using the light for a week but haven't noticed any benefit, please call your Lumie distributor or speak to your doctor.

IMPORTANT: the guidelines in this leaflet are not medical directions for treating any condition. If your doctor has suggested using light therapy, follow their instructions and discuss any effects with them.

3 Safety

- Keep the unit away from water and damp. For indoor use only.
- Use only the Lumie adaptor supplied.
- Children using the light, or playing near it, should be supervised.
- Do not attempt to service any part of the unit.
- This product conforms to the electromagnetic compatibility standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens simply move the product away from the source of interference.
- At the end of its life, do not throw the unit out with normal domestic waste. Please take it to an official recycling point to be disposed of responsibly.

4 Getting started

If the product is damaged please contact us before use.

Position Desklamp at arm's length (about 50cm/20in) with the head angled towards your face. Connect the power adaptor and plug into the mains. You can turn the light on and off and adjust the brightness using the touch-sensitive area on the base.

Touch control lighting

The 'sense' touchpad provides adjustable lighting that means you can use Desklamp for reading and task lighting as well as light therapy.

There are four light levels offering approx. 40, 60, 80 and 100% of full brightness. Tap the 'sense' touchpad to turn on the light at low level, to increase the brightness and then to turn off the light.

Or apply a continuous touch to cycle through the four light levels and off - simply let go at the desired light level.


5 Using Desklamp

Desklamp has a removable diffuser that easily lifts off to provide flexible light therapy as well as excellent task lighting.

Position Desklamp at arm's length (about 50cm/20in) with the head angled towards your face. You shouldn't stare at the light – especially when using Desklamp without the diffuser – but light must reach your eyes to have any effect. If you find your eyes ache slightly, try using it further away, or with the light dimmed or switch it off for a while.

Without the diffuser – typically 30 minutes

With the diffuser off and the light at full brightness, you're getting the maximum benefit from Desklamp so you don't need to use it for long. This provides bright, intensive light therapy (2,500lux at 50cm) that's ideal if you're pushed for time or need a quick boost at work.

With the diffuser – typically 60 minutes

With the diffuser in place and the light at full brightness, Desklamp provides a softer, wider spread of light (1,250lux at 50cm) that's more suitable if you have time for a longer session.

There's no harm in using Desklamp for longer than 60 minutes, e.g. if you're away from your desk you can make up this time by dipping in throughout the day. If you're more comfortable with it further away or with the light at less than full brightness then you'll need to use it for longer.

You can use Desklamp whilst you're working at the computer, watching TV, studying, crafting etc. Dark glasses, tinted lenses or closing your eyes will all limit the effects of light therapy.

You should notice positive benefits after using Desklamp for three or four days in a row. You can monitor your progress using the questionnaire at the end of the instructions.

6 When should I use it?

It's best to tackle the winter blues or SAD as soon as possible. Look out for early signs - often in September or October - and if you start to feel lethargic, anxious or irritable, start using your Desklamp.

Most people can skip a day here and there, especially when the weather is bright. As spring comes around you will probably start to cut back your light therapy; your body will tell you if it's too soon and if you notice your symptoms returning simply use your Desklamp again for a few days.

Light therapy is proven to put you in a better mood, to boost alertness and to help you feel more awake so Desklamp is ideal for use at work. Use it wherever and whenever you need to restore concentration and give yourself a mental boost.

7 What time of day?

If you're using Desklamp to help with the winter blues or SAD, bright light in the morning has been shown to be most effective so try that first. For some people, early evening light exposure is beneficial. Whether you're using Desklamp for winter blues, SAD or as a general pick-me-up, once you've found a pattern that's right for you, try to use it around the same time every day.

I find it difficult to get going in the morning

Switch on as soon as you can. If you don't have much time in the morning, then you should aim for a top-up or full session in the afternoon. You could take your Desklamp to work if you're at a desk every day.

I feel sleepy or fall asleep too early

Use Desklamp in the late afternoon/early evening. It's best to avoid bright light too close to bedtime - within about three hours - as you may find it difficult to get to sleep.

8 LED lifetime

LEDs are different from light bulbs in that (i) they very gradually decay rather than suddenly stop emitting light and (ii) they cannot be simply replaced by another LED unit. LED manufacturers usually express the lifetime of LEDs as the time to reach 70% of their initial light output and this figure typically exceeds 50,000 hours. This is at least five years and can be 20 years depending on usage. However, it is important to also bear in mind that LEDs can spontaneously fail.

9 Warranty

This unit is covered under warranty against any manufacturing defect for three years from the date of purchase. This is in addition to your statutory rights. If you are returning your Desklamp for service, it must be in its original packaging, if possible, including internal pieces. Make sure that the plug is prevented from moving in transit and damaging the unit. If the unit is received damaged, then we may have to charge to repair it. Please call Lumie for further advice.

10 Technical specification

Subject to change without notice.

Lumie Desklamp

Lumie 100-240Vac, 50/60Hz power adaptor (supplied) 

Operating temperature: 0°C-35°C

Class IIa Medical Device

11 Contact us

Lumie is Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits. If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on +44 (0)1954 780500 or email info@lumie.com.

Lumie
3 The Links, Trafalgar Way
Bar Hill
Cambridge
CB23 8UD
UK

Lumie is a registered trademark of:

Outside In (Cambridge) Limited, incorporated in England and Wales.

Registered Number: 2647359. VAT No: GB 880 9837 71.

© Lumie 2012

This questionnaire is designed to help you assess how much the lights are helping. It is based on the standard tool used by doctors. Not everyone will have every symptom but your total score should go down over the weeks. You might like to talk it over with your doctor or another person in your household.

Compared with when you are well, how have you felt during the past week?

0 = no different 1 = a little bit 2 = somewhat 3 = quite a lot 4 = definitely, badly

DEPRESSION SYMPTOMS CHART					
I have been...	before using the light	after 1 week	after 2 weeks	after 3 weeks	after 4 weeks
down and depressed					
less interested in doing things					
less interested in sex					
less interested in eating					
losing weight					
difficulty getting to sleep					
sleeping badly, waking at night					
waking too early					
feeling a failure					
feeling that life is not worth living					
tense, irritable, worried about things					
convinced I'm ill, excessively worried about my health					
thinking or speaking more slowly					
fidgety, restless, can't keep still					
feeling cut off from what's going on, as if in a dream					
paranoid, suspicious					
preoccupied, easily distracted					
physically tired e.g. heavy limbs, headaches, aching back/muscles					
other physical symptoms e.g. stomach problems, sweating, urinating frequently, dry mouth, cramp, heart palpitations, hyperventilating					
morning is worse than evening					
evening is worse than morning					
Standard depression score					
feeling unsociable, avoiding people					
gaining weight					
wanting to eat more than usual					
ACTUALLY eating more than usual					
craving sweet and starchy foods					
sleeping too much, feeling more tired than I should be					
worse in the afternoon or evening, then better at least an hour before bedtime					
Supplemental depression score					
TOTAL standard and supplemental score					

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Outside In (Cambridge) Ltd